

# Biological indicator OPS for disinfection processes of surgical clogs

## Product information

**Field of application:** OPS is a biological indicator, which is designed for the validation and the routine control of cleaning and disinfection processes of washer disinfectors for medical clogs. (> 60 °C)

**Features:** OPS indicators contain populations of *Enterococcus faecium* and are contaminated with test soil according to ISO 15883-7. OPS indicators are suitable for the qualitative verification of the required log 5 reduction.

**Specifications:** *Organism: **Enterococcus faecium***  
*Mean population (cfu): **≥ 10<sup>5</sup>***  
*Carrier material: **stainless steel (appr. 140 x 10 mm)***  
*Primary packaging: **paper / foil***  
*Organic burden: **RAMS and sheep blood***  
*Shelf life: **3 months from the date of manufacturing***

**Storage:** + 4 °C to + 8 °C

**Disposal:** After disinfection process, dispose with domestic waste.

**Packing Unit:** 10 pcs

**Order No:** BI-OPS-19001

### Example of use:

1. Take the biological indicators (carrier with test soil) out of the pouch and attach them with cable ties firmly to the medical clogs or if this is not possible to the collecting basket. The contaminated side has to face outside. One indicator is meant to be a growth and transport control. Do not wash the control indicator.
2. Sanitize your hands, once the indicators are fixed.
3. Bring the wash rack into the washer disinfectant. If possible turn off the drying program because hot air drying may falsify the result. Check selected program and start the program.
4. When the disinfection process is finished, take the indicators aseptically out of the WD. In case there are no sterile tweezers available use sterile one-way gloves. Cut the cable ties. Make sure you touch the indicators only on the outer edge close to the drill hole.
5. **Put each indicator in a separate tyvek-pouch and seal it. Sanitize your hands before every transfer of the next indicator.**
6. Incubation: 4 days at 35 °C ± 2 °C (e.g. incubate with an *Enterococcus* selective nutrient media)
7. Daily check all tubes for growth of the test organism.
8. Note down the results. The results are only valid if the growth control shows typical growth.