

PurePep® EasyClean Purification Kits

GYROS PROTEIN
Technologies

Starter Kit

Our Starter Kits are the perfect introduction to the world of PurePep EasyClean (PEC) peptide purification technology. The ready-to-use package comes with PEC-Linker, cartridges prefilled with activated filter material, and related consumables. You can choose between two formats that provide you with the ideal kit in a suitable number and scale for your purification needs.

PurePep EasyClean Starter Kit 8 x 25 µmol*
meet your research demands

PurePep EasyClean Starter Kit 8 x 100 µmol*
tackle the elevated scales

*referring to synthesis scale



“ The time saved working with the Starter Kit is precious, enabling technicians to work on other tasks. Also, we were offered fantastic customer support, with fast communication and great suggestions on how to solve issues to achieve the highest yield and purity.”


Senior Scientist, Pharma Company

READY-TO-USE



Product Highlights:

- The simplest access to PEC-grade peptides - manual or automated setup
- Fastest access to **8 or more purified peptides** per day
- **25 to 100 µmol** R&D Scales – **up to 200 mg** of product per syringe
- **75% lower** solvent consumption and waste generation
- Reduced purification cost



PEC-Linker RC+

	PurePep EasyClean Starter Kit 8 x 25 µmol	PurePep EasyClean Starter Kit 8 x 100 µmol
Scale	8 x 25 µmol	8 x 100 µmol
PEC-Linker	RC+	RC+
Activated filter material	Agarose100 filled in 8 fritted cartridges	Agarose100 filled in 8 fritted cartridges
Reducing Agent	1,4-Dithiothreitol	1,4-Dithiothreitol
Blocking Agent	L-Cysteine	L-Cysteine
Buffer	Mixture of citric acid-sodium carbonate	Mixture of citric acid-sodium carbonate
Optional equipment	Vacuum manifold "Starter" (P0020920)	Vacuum manifold "Starter" (P0020920)
Documentation	Manual in English	Manual in English
Product code	P0020873	P0020872