

90XL Meter Ordering Guide

Contact us: mesalabs.com/renal-care



● Feature Available ○ Feature Unavailable

Products	System 5521106000	Meter-CT/ Press 5521106001	Kit 5521106002	Meter-2 CT 5521106007
90XL Display Module	●	●	●	●
90XL User Guide	●	●	●	●
Cover/Holder with IV Pole Clamp	●	●	●	●
Conductivity/Temperature Module with integrated cable	●	●	●	2
Pressure Module with integrated cable	●	●	●	○
pH Module with integrated cable	●	○	○	○
Negative Pressure Adapter	●	○	●	○
In-line Adapter	●	○	●	○
Transducer Protectors	2	2	2	2
Monitoring Line with stopcock	2	●	●	●
90XL 110V Adapter	●	●	●	●
90XL 110V Blade	●	●	●	●
Sampling Kit	●	○	○	○
Standard Solutions (see below)	●	○	○	○
Carrying Case	●	○	○	○
Standard Solutions				
1.0 mS Conductivity Calibrator Solution, 16oz./475ml	●	○	○	○
14.0 mS Conductivity Calibrator Solution, 16oz./475m	●	○	○	○
100 mS Conductivity Calibrator Solution, 16oz./475m	●	○	○	○
4.0 pH Buffer Solution, 16oz./475m	●	○	○	○
7.0 pH Buffer Solution, 16oz./475m	●	○	○	○

Payment Terms: NET 30 WAC, MasterCard, VISA and AMEX Orders Accepted

Shipping: Prepaid and added to Invoice with a nominal handling charge. *Product specifications and pricing are subject to change without notice.*

Rev. 12.2022