

Universal Torqo Chuck

Torque analyzers are customized for the specific bottle and cap combination to be used on the unit. Chucks are most often designed to match the serration patterns of each cap resulting in additional change parts, complicated changeover processes and greatly increased cost to the end user.

As an increased variety of closure designs have gained prominence and popularity in the industry, Mesa Labs has developed a new Universal Chuck to accommodate multiple serration patterns and cap sizes with a single change part.



Varying serration patterns for similar cap diameters

Features

- Fits any Torqo model 1502 – 1600 torque testers
- ROPP cap compatible
- Tapered fit to accommodate most 26–32mm OD caps (CR, CT, TEB, etc)*
- Compatible with various serration patterns
- Eliminates the need for multiple chucks
- Non-destructive (does not leave marks when used properly)
- Compatible with smooth caps (torque limitations apply)

Contact us today for pricing info or to inquire about universal chuck designs for other cap OD ranges.

**The composition of closure, serration pattern, upper torque specification limit and bottle shape may influence compatibility. Before placing an order, please send sample closures to our customer service for verification.*