

Product Specification

Population Assay Materials

Product Description

Population Assay Materials include glassware and growth media, which are sold separately. Glassware (PAK-G) includes tubes, pipettes, and beads required to make appropriate dilutions following the Mesa Labs Population Assay procedure. The media is sterile Difco™ brand Tryptic Soy Agar appropriate for growth of spores. This is the same media used by Mesa Labs for population assays.

Technical Parameters

Media	Shelf-life	12 months from date of manufacture.
	Storage	Room Temperature
	Volume	240 mL
Glassware	Tubes	x4 19.5 x 145mm flat bottom with four 6mm glass beads each
		x12 16 x 125mm dilution blank
	Pipettes	x2 10mL
		x2 5mL
		x8 2mL
		x8 1mL

Indications for Use

Mesa's Population Assay Materials contain everything required to perform accurate population assays of Mesa Labs' biological indicator products.

Packaging


PAK-G glassware is packaged in a corrugated box. The PAK-M media is available in bottles of 240mL.

Labeling

- Each media bottle is labeled with lot number and expiration date, with additional information on the box label.
- Each PAK-G box has a self-adhesive label on the outside to identify the product.

Reorder Information

- PAK-G (Glassware and pipettes)
- PAK-M (Difco Growth Media)



PAK-M



PAK-G