



## PRODUCT CHANGE NOTIFICATION

29 Nov 2018

**Notification ID:** CN-181103

**Notification Type:** Product Specification Change

Dear Customer,

We would like to take this opportunity to notify you of a specification change to Mesa Labs **ACD/6 MesaStrip Biological Indicator for Chlorine Dioxide Sterilization**. Mesa Labs no longer supports a 48-hour reduced incubation time (RIT) claim for ACD/6 when cultured in RM/100 Mesa Labs Red Releasat Culture Media.

**Please note this change does not affect use of this product when incubated for 7 days.**

The change was made because Mesa no longer has a chlorine dioxide test vessel which is necessary to continue to support this RIT label claim.

As always, if you have any questions, please contact your Mesa Laboratories Representative.

  
Beth Ridgeway, Director – Manufacturing

24 Nov 2018  
Date

  
Janis E. Smoke, Director – Quality & Regulatory Affairs

29 Nov 2018  
Date